

Non-Autoclavable Adjustable Pipette



Autoclavable Adjustable Pipette



Bottle-Top Dispenser



Motorized Pipette Filler

Liquid Handling Devices

Specifications and equipment are subject to change without any notice or obligation on the part of the manufacturer. **Copyright © CATALYST BIOTECH**. PLEASE READ THE CORRESPONDING MANUALS CAREFULLY BEFORE USING YOUR EQUIPMENT.



Catalysi Biotisch

C-302, Third Floor, Anupam CHS Ltd. Plot No.: 79, Old Panvel - 410 206.

Maharashtra - India.

E-mail: catalyst@catalystbiotech.com Website: www.catalystbiotech.com Our Authorised Distributor :





CatPette Bottle-Top Dispensers are fully autoclavable at 121°C and are ideal for safe, accurate and precise dispensing of bases, saline solutions, acids and some organic solvents. It's equipped with easy to clean and maintain tools.



Smooth piston design provides effortless dispensing

Quick and precise volume adjustment

Chemical resistant to most aggressive liquids and acids

Free of additional adapters, allows fitting to standard reagent bottles.

Anti-drip cap prevents accidental discharge.

Closure cap protects, against reagent contact and is easy to attach and remove, even while wearing gloves

Vapour pressure	Max. 500mbar		
Viscosity	Max. 500mm²/s		
Temperature	Max. 40°C		
Density	Max. 2.2g/cm³		



CatPette

Four models are available, covering volumes between 0.5ml to 50ml and they come with 5 adapters (45/40mm. 45/38mm, 45/32mm, 32/28mm & 32/25mm) to fit all common lab bottles.



CatPette: Bottle-Top Dispenser							
Cat. No.	Volume Range	Increment	Accuracy ≤ ±		Coefficient of Variation ≤		
			%	μΙ	%	μl	
1030101000	0.5 - 5ml	0.1ml	0.5	25	0.1	5	
1030102000	1.0 - 10ml	0.2ml	0.5	50	0.1	10	
1030103000	2.5 - 25ml	0.5ml	0.5	125	0.1	25	
1030104000	5.0 - 50ml	1.0ml	0.5	250	0.1	50	

